


# **Dental Implant Brand Detector**

An Al-based solution that makes dentists work more efficiently





## **Meet Our Team**



**Advisor** 





Zongyuan Yu




Jack Chen




Clinical Assistant Professor

Utku Pamuksuz








Zoey Zhou



Fiona Fei



# **Executive Summary**

- We developed an Al solution which takes a panoramic x-ray and identifies the implant brand for dentists in real time.
- Though the data is limited, our Al detector achieves a desirable accuracy.
- Our Al ensemble solution is estimated to save 30 minutes per consultation session for dentists.
- The App is ready to deploy!







# A Dental Implant is a replacement of tooth root



- The dental implant surgery is a procedure that replaces tooth roots with metal, screw-like posts and replaces damaged or missing teeth with artificial teeth that look and function much like real ones.
- The average lifespan of a dental implant is approximately 10 to 15 years

#### Structure of The Dental Implant System







## **Business Problems**



❖ For patients who have implants already and require further treatment or replacement, identifying brands of the existing implants on images is a vital task to the dentists

#### **Problems**

- Patients usually don't know the brand of their dental implants.
- Even if the patients know their brands, they are likely to forget the name of brand years later when a replacement is needed.





# Deep Learning based approach could replace the manual searching process



#### **Current solution**

Dentists are currently searching for the origin of the implant manually and

developing a treatment plan once they find it.

#### **Desired solution**

- Automate this process to save medical professionals' time
- Solution should automatically detect implants in an x-ray image and predict the brands of the implants





# Data are X-Ray Images including Panoramic and Bitewing

**Providing Company:** Great Lake Dental




Data Format: X-ray Images, YOLO labeled data


#### **Examples:**

Panoramic




Bitewing







# 8 Classes are Labeled and 'Existing-Implant' is the main focus



#### YOLO (You Only Look Once) Labeled Data:

#### 8 Classes in Total:

Class Name	Label Number	
existing - filling	0	
existing - root canal treatment	1	
existing - implant	2	
caries low risk	3	
caries moderate risk	4	
caries high risk	5	
periodontitis	6	
crown	7	

#### **Ready for Processing:**

Almost 2 million panoramic X-Rays/bitewing/fmx/photo

### **Annotated & Ready for Training:**

30 thousand panoramic X-Ray

#### **Scheduling Optimization:**

400 implant labeled images along with Dental Clinic's database records with 11 brands in total







# **Solutions to Three Main Challenges on Data**

1. Given limited number of labeled implants images (<400) available, it is not sufficient for deep learning algorithms.

Solution: Use YOLO to label the unlabeled images to get more samples

A lot of X-ray images are blurry or overexposed.
 Solution: Manually throw out bad quality images

 Dentists will reshape each implant before planting. So it will look different even though they might be exactly the same implant.
 Solution: Manually categorize the implants




Figure 1.Blurry, Overexposed implant images





Figure 2. Hard to tell whether same or different brand



# The Analytical Plan Follows a Stepwise Manner and Attempt Different Techniques

Model for detect and label the implants locations

Build the final clustering model with the complete dataset





Compare different models and evaluate their performance







# **YOLO Labeling technique**









#### **Procedures:**

- 1. Use YOLO to detect each tooth in the X-ray
- 2. Iterate through every tooth detected
- 3. Label each tooth using a box for future analysis





# The size of available image data is tripled by implementing YOLO for implant detection

To filter the data and locate implants, we train YOLO model with known implant images. Then use the trained model to find the unlabeled implants.

Use **Intersection Over Union (IOU)** as our metric to train the model, which is Area of overlap between predicted bounding box and true bounding box

Increased our available implants image data from 390 to 996.





Figure 1. X-ray Images Before Labeling





Figure 2. X-ray Images After Labeling



# **11** Classes are Manually Categorized for the Training Purpose





























We manually categorize our implants into 11 brands, trying to gather as many data as possible

We have a total of 933 implant images in 11 classes. However, the data are imbalanced






# Models 1: K-Means Clustering did not generate promising results

#### **#1 K-Means Clustering**

#### Procedure:

- 1. Feature extraction with pretrained VGG16
- 2. Use K-Means clustering model with 11 classes.



-20 --30 -20 -10 0 10 20 30

Figure 1. Cluster colored with true classes

Figure 2. Cluster colored with predicted classes

#### Results:

It did NOT give us a good cluster.

The cluster on the right (Cluster from K-Means) is very different from the cluster on the left (True Class)

**Decision:** 

NOT to go forward with it



#### Models 2: EfficientNet was associated with low accuracy and imbalanced issue

#### #2 Pre-trained EfficientNet Model

#### Procedure:

Use Pre-trained EfficientNet to predict classes
Apply data augmentation that create more sample by flipping, resizing the images.

#### Results: 23% Accuracy with heavy classes imbalance

On the left (Predicted Results), only class 4 and 9 have lots of correct prediction which is due to over 50% of our training samples are from class 4 and 9.




Figure 1. Confusion Matrix — Most of the predictions are class 4 and 9

Classification Report					
	precision		recall	f1-score	support
	01	0.00	0.00	0.00	9
	02	0.00	0.00	0.00	16
	03	0.00	0.00	0.00	8_
	04	0.49	0.32	0.39	84
	05	9.99	<b>0.00</b>	9.99	8
	96	0.00	0.00	0.00	6
	97	0.00	0.00	0.00	4
	08	0.00	0.00	0.00	5
	09	0.17	0.53	0.26	30
'	10	0.00	0.00	0.00	14
	11	0.00	0.00	0.00	3
2001	2261			0.23	187
accur		0.06	0.00		
macro	_	0.06	0.08	0.06	187
weighted	avg	0.25	0.23	0.22	187

Figure 2. Classification Report — Class 4 and 9 has much better precision and recall than other classes



# Models 3 : Siamese Network is more robust to class imbalance and require less training data

#### **#3 Siamese Network**

Siamese Network is introduced since it is able to use very small number (>6) of images to achieve accurate predictions

#### Model Feature:

Siamese Network, instead of learning all the features of the classes, learns to differentiate between classes




Figure 1. Siamese Network - Signature Verification Example





#### **Models 3: Dissimilarity Scores from Siamese Network**

#### **#3.1 Siamese Network example output**

Siamese Network only provides similarity score between two images.


#### Following steps to predictions


- 1. Compare a test image to every train images we have and get the dissimilarity scores for all such pairs.
- 2. Get an average dissimilarity score grouped by classes.
- 3. The class with the best score will be the top prediction.





## **7** Classes are Manually Categorized for the Training Purpose





















#### Distribution of Classes of Implant



Every implant from different classes have distinct characteristics.
Untrained human can tell apart implants that are from different classes.

Total 143 images gives us ~10000 unique image pairs for training Siamese Network






### Models 3: Re-apply Siamese Network with 7 Manually Categorized Classes

#### **#3.2 Siamese Network Evaluation**

#### Recap on how we make prediction:

- Compare a test image to every train images we have and get the similarity scores for all.
- 2. Get an average dissimilarity score from each class.
- 3. Classes with the best scores will be the top predictions.

Classes	1	2	3	4	5	6	7	Overall
Accuracy	47%	5%	5%	9%	72%	38%	14%	31.5%
Top 2 Accuracy	77%	23%	5%	9%	91%	95%	14%	58.7%



#### **Metrics and evaluations**

#### **Metrics:**

- 1. For YOLO model, Intersection Over Union (IOU), which is a standard metric for evaluating object detection algorithm, is used to evaluate our YOLO model
- 2. For multiclass classification model, we use both accuracy and top-2 accuracy to evaluate our clustering model


Metrics	
Intersection Over Union (IOU)	Area of overlap/Area of union (between predicted bounding box and true bounding box)
Accuracy	Out of all predicted as class A, how many really belong to class A?
Top-2 Accuracy	If the true brand of the implants is within our predicted top 2 most probable brands, we consider it as correct prediction.







# **Final Product Demo**





# **Final Conclusion**



#### Conclusion

Our data was not enough for training the K-Means or pre-trained deep learning model. Siamese Network could overcome this problem and performed the best out of all models.

#### **Deliverable**

Fully functional implant identification app based on Siamese Network model

#### **Next Step**

- In the future, the model's accuracy could be improved if we have access to higher quality data (implants' brands labeled by medical professionals, etc.)
- 2. Similar model could be used for training other Al for dentists. (e.g., detecting the existing tooth filling material when operating a tooth filling replacement)





#### Contrastive Loss Function where

- Dw is is defined as the euclidean distance between the outputs of the sister siamese networks
  - Y indicating whether the two inputs are from the same class (0 means same classes)
- m is a margin value which is greater than 0. Having a margin indicates that dissimilar pairs that are beyond this margin will not contribute to the loss:

$$(1-Y)\frac{1}{2}(D_W)^2 + (Y)\frac{1}{2}\{max(0, m-D_W)\}^2$$